

VENOSAN® Compression Stockings

Product guide



Simple. Effective. Swiss quality.

VENOSAN® Compression Stockings help to make life easier in many situations. The pressure VENOSAN Compression Stockings exert promotes healthy veins. They can help to prevent venous disease as well as being an important part of the treatment for patients being treated for acute and chronic venous disorders.

The diverse VENOSAN range offers the right product for a wide variety of needs. All VENOSAN Compression Stockings have the following characteristics in common:

Simple.

The VENOSAN range provides a simple guide to help you determine the correct size, compression class and style of compression stocking.

The wear and care instructions included with VENOSAN Compression Stockings are easy to follow which patients appreciate.

Effective.

VENOSAN Compression Stockings are RAL-certified, meaning that they provide effective therapeutic compression. The effective design of VENOSAN Compression Stockings means that for patients, many will experience enhanced treatment success and improved quality of life.

Swiss quality.

Products from Switzerland are known for excellent workmanship and top quality.

Manufactured using high-quality materials, such as the original Lycra® and elastane, VENOSAN Compression Stockings offer wearing comfort and durability.



Product overview

Medical compression stockings and pantyhose.

VENOSAN® 4000

CLIMATE EFFECT

- Climate effect
- Semi-transparent
- Comfortable
- Highly elastic


Indications

- CVI (chronic venous insufficiency) 1
- Varicosis during pregnancy
- After vein surgery



Composition

71% TACTEL® polyamide
29% elastane



VENOSAN® 5000

THE MAGIC OF SILVER

- Odour-inhibiting
- Anti-static
- Temperature-regulation
- Anti-bacterial

Indications

- CVI (chronic venous insufficiency) 1–2
- After vein surgery
- Mild skin conditions
- Lipoeidema stage 1



Composition

40% TACTEL® polyamide
30% elastane
20% cotton
10% polyamide/silver (X-Static®)



VENOSAN® 6000

CLASSIC

- Maximum sheerness
- Breathable
- Ease of care


Indications

- CVI (chronic venous insufficiency) 1-2
- Varicose veins
- Varicosis during pregnancy
- After vein surgery



Composition

70% nylon polyamide
30% elastane



Product overview

Medical compression stockings and pantyhose.

VENOSAN® 7000 NATURAL TOUCH

- High stiffness (resistance)
- Skin-friendly
- Temperature regulation

Indications

- CVI (chronic venous insufficiency) 2–3
- Lipoedema stage 2–3
- After vein surgery
- After thrombosis
- PTS (post-thrombotic syndrome)
- Mild lymphoedema
- Marked orthostatic stress



Composition

62% TACTEL® polyamide
26% elastane
8% cotton
4% Lyocell/SeaCell®



VENOSAN® 8000 ULCERFIT SUCCESSFUL THERAPY

- Two-layer stocking system
- Fast and easy to use
- High wearing comfort
- Liner can be worn at night

Indications

- Treatment of venous leg ulcers
- Helps to prevent recurrence of venous leg ulcers



Composition

Overstocking:
62% TACTEL® Polyamide
26% Elastane
8% Cotton
4% Lyocell/SeaCell®

Liner:

90% Polyamide
10% Elastane



VENOSAN® SOFT VENOSAN® PRO VENOSAN® STRONG FLAT KNIT THERAPY

- Flat-knit custom made garments
- Smooth and durable
- Low to high stiffness
- Manufactured to European Standards

Indications

- Treatment and management of Chronic Oedema and Lymphoedema



Composition

Polyamide & Elastane









VENOSAN PRO (CCL II)
& VENOSAN STRONG (CCL 4S)
also contains cotton

VENOSAN PRO also contains Latex







Product lines

Medical compression stockings and pantyhose.

										
Line	Colours	CCL	AD Below knee stocking	ADH Below knee stocking with dotted top band	AG Thigh- length stocking	AGH Thigh-length stocking with dotted top band	AGH Thigh-length stocking with lace top band	AGG Thigh- length stocking with belt	AT Pantyhose	ATU Maternity pantyhose
VENOSAN® 4000 CLIMATE EFFECT	Mexico	CCL 1	●	●	○	●	●	○	●	●
	Marokko	CCL 2	●	●	○	●	●	●	●	●
	Black	CCL 2	●	●	○	●	●	●	●	●
VENOSAN® 5000 THE MAGIC OF SILVER	Mexico	CCL 1	●	●	○	●	●	●	●	○
	Marokko	CCL 2	●	●	○	●	●	●	●	○
	Black	CCL 2	●	●	○	●	●	●	●	○
VENOSAN® 6000 CLASSIC	Beige	CCL 1	⊙			⊙			⊙	
	Black	CCL 2	⊙			⊙		⊙	⊙	
		CCL 3	●	○	○	○	○	○	○	○
VENOSAN® 7000 NATURAL TOUCH	Beige	CCL 1	⊙							
	Black	CCL 2	●	●	○	●	○	●	○	
		CCL 3	●	●	○	●	○	○	○	
VENOSAN® 8000 ULCERFIT SUCCESSFUL THERAPY	Beige	CCL 2	⊙							
		CCL 3	⊙							
VENOSAN® SOFT VENOSAN® PRO VENOSAN® STRONG FLAT KNIT THERAPY Contact your local distributor for a full list of available garments and compression classes	Beige	CCL 1	○	○	○	○		○	○	○
	Black	CCL 2	○	○	○	○		○	○	○
		CCL 3	○	○	○	○		○	○	○
		CCL 4	○	○	○	○		○	○	○

 available in standard sizes
  available in custom sizes only
  available in standard and custom sizes

Pressure at ankle  18–21mmHg
  23–32mmHg
  34–46mmHg
  >49mmHg

Accessories for compression stockings.



VENO-GLIDER

Foot Slip Donning Aid

Article no.	Description	Unit
995 345	VENO-GLIDER	1 pc.



VENOSAN® Dot Gloves

Washable gloves for donning and doffing compression stockings, available in two sizes.

Material composition

Knitted fabric: 100% nylon polyamide,
Dots: PVC

Article no.	Description	Unit
VC 003	VENOSAN Dot Gloves Size S/M	1 pair
VC 004	VENOSAN Dot Gloves Size L/XL	1 pair



Socks-jet

Donning and doffing aid
for compression stockings

Article no.	Description	Unit
AV 002	Socks-jet with handles	1 pc.



COLTEX

Special adhesive for securing medical
compression stockings to the skin.

Composition

PVP/VA copolymer, glycerin,
2-propanol

Article no.	Description	Unit
KC 67000	COLTEX Adhesive 60ml	1 pc.

For more information, please contact:

Complete Medical Compression WA

P: 08 6174 0314

M: 0425124 576

E: sales@cmompression.com.au