### **ExoStrong**<sup>™</sup> Affordable alternative to custom garments.

ExoStrong is a full line of off-the-shelf flat knit gradient compression garments for the maintenance of lymphoedema or oedema during the day or periods of high activity.

For lymphoedema patients, flat knit garments are one of the primary self-care options. Where circular knit migrates into the skin folds, flat knit helps shape the limb and provides the high containment needed to maintain reductions made in therapy.

ExoStrong is a fraction of the price of custom flat knits allowing patients to choose the right daytime garment at an affordable price.



"ExoStrong is a great option for my patients that fit into standard sizes, because it provides them with the containment and comfort of flat-knit — at an affordable price compared to custom made."



#### Made in the USA

High-quality flat knit, manufactured at LRUS in Milwaukee, Wisconsin.



#### **Off-the-Shelf**

Complete range of garments at a lower cost than custom-made options.



#### Supports mobility

High containment garment in a comfortable and functional design.



#### Quick and easy fitting

Only 2–5 measurements are needed to select your size for each style.



#### Easy garment care

Machine washable and dryable to maintain fit for the life of the garment.



#### Service guarantee

Guaranteed to be free of material and manufacturing defects for the six-month life of the garment.

## Lohmann & Rauscher

#### **ExoStrong**<sup>™</sup>

#### Off-the-shelf flat knit compression garments.

#### Indications

- Mild to severe lymphoedema
- Oedema

#### Intended Uses

- Helps maintain reduced volume and shape by supporting extremities
- Management of mild to severe lymphoedema
- Management of mild to severe oedema (non-blood related)

#### Contraindications

- Arterial insufficiency/degeneration or ABPI < 0.8</p>
- Moderate to severe peripheral arterial disease
- Untreated congestive heart failure
- Untreated/systemic infection (i.e. cellulitis)
- Allergies or sensitivities to materials

#### Precautions

- Impaired sensation (peripheral neuropathy)
- Cognitive impairment

#### Care Instructions

#### Washing

- Machine wash warm on permanent press
- Do not use bleach

#### 

Do not dry clean

#### Material Content

Nylon (polyamide), lyocell, spandex (elastane), silicone\*

\*ExoStrong garments with Silicone Top only. Latex free.

For all queries relating to Exo products as well as ReadyWrap<sup>®</sup>, SwellSpots<sup>™</sup>, Caresia<sup>™</sup>, TributeNight<sup>™</sup>, Tribute<sup>®</sup> Wrap, or ExoSoft<sup>™</sup> products please contact: **Complete Medical Compression WA Phone:** 08 6174 0314 Email: sales@cmcompression.com.au


Caresia™, ExoStrong™, ExoSoft™, SwellSpots™ and Tribute Night™ are trademarks of Lohmann & Rauscher GmbH & Co KG. © 2024 ReadyWrap® and Tribute® are registered trademarks of Lohmann & Rauscher GmbH & Co KG. © 2024

Lohmann & Rauscher Pty Ltd · 321 Ferntree Gully Road · Mount Waverley, Victoria 3149, Australia · Customer Service: 1300 572 869 Lohmann & Rauscher Limited · 21 Queen Street · Auckland 1010 · New Zealand · Customer Service: 0800 572 869





## **ExoStrong**<sup>™</sup> Off-the-shelf flat knit compression garments.



Tumble drv low Do not iron

Drying

Scan QR Code to view Upper extremity

instructions video

## **ExoStrong**<sup>™</sup> Upper Extremity garments



# Gauntlet

sold individually Compression 20-30mmHg

	Circumfer	ence (cm)	REF				
	Palm (cB)	Wrist (cC)	Beige	Black			
Small	17 – 19	15 – 17	104296	104297			
Medium	19 – 21	17 – 19	104298	104299			
Large	22 – 24	20 – 22	104300	104301			
X-Large	25 – 27	23 – 25	104302	104303			



#### Glove

sold individually Compression 20-30mmHg

					Finger	Length				
Size	Circumfe	rence (cm)	Quarter		Half		Full			
3120	Palm (cB)	Wrist (cC)	Beige	Black	Beige	Black	Beige	Black		
	ram (cb)		REF							
Small	17 – 19	15 – 17	104272	104275	104273	104276	104274	104277		
Medium	19 – 21	17 – 19	104278	104281	104279	104282	104280	104283		
Large	22 – 24	20 – 22	104284	104287	104285	104288	104286	104289		
X-Large	25 – 27	23 – 25	104290	104293	104291	104294	104292	104295		

#### Armsleeve

sold individually Compression 20-30mmHg

		Top of Armsleeve						
Size -	Length (cm)	Circumference (cm)			Kı	nit	Silicone	
	(10, 0)	Wrist	Elbow	Axila (cG)	Beige	Black	Beige	Black
	(IC-G)	(cC)	(cE)		REF			
Small	Short (up to 43cm)		24 – 27	27 – 30	104304	104305	104306	104307
	Average (up to 46cm)	15 – 17			104320	104322	104321	104323
	Tall (up to 49cm)				104336	104338	104337	104339
Medium	Short (up to 43cm)		27 – 31	31 – 35	104308	104310	104309	104311
	Average (up to 46cm)	17 – 19			104324	104326	104325	104327
	Tall (up to 49cm)				104340	104342	104341	104343
	Short (up to 43cm)	20 – 22	31 – 35	35 – 40	104312	104314	104313	104315
Large	Average (up to 46cm)				104328	104330	104329	104331
	Tall (up to 49cm)				104344	104346	104345	104347
- X-Large	Short (up to 43cm)		36 – 41	40 – 45	104316	104318	104317	104319
	Average (up to 46cm)	23 – 25			104332	104334	104333	104335
	Tall (up to 49cm)				104348	104350	104349	104351

All measurements in centimetres.

## **ExoStrong**<sup>™</sup> Lower Extremity garments

#### **Below Knee**

sold individually Compression 20-30mmHg Top Style: Open

							Top of Stocking			
Size	Length (cm)			Circumference (cm)		Knit		Silicone		
Size		(1.6.)	(ID)	(cB)	(cC)	Beige	Black	Beige	Black	
		(IA)				REF				
	Short	up to 17	up to 38		27 – 36	104181	104182	104183	104184	
Small	Average	up to 19	up to 41	18 – 20		104201	104203	104202	104204	
	Tall	up to 22	up to 44			104221	104223	104222	104224	
	Short	up to 17	up to 38	20 – 24	29 – 39	104185	104187	104186	104188	
Medium	Average	up to 19	up to 41			104205	104207	104206	104208	
	Tall	up to 22	up to 44			104225	104227	104226	104228	
	Short	up to 17	up to 38	24 – 28	32 – 43	104189	104191	104190	104192	
Large	Average	up to 19	up to 41			104209	104211	104210	104212	
	Tall	up to 22	up to 44			104229	104231	104230	104232	
	Short	up to 17	up to 38		33 – 45	104193	104195	104194	104196	
X-Large	Average	up to 19	up to 41	28 – 31		104213	104215	104214	104216	
	Tall	up to 22	up to 44			104233	104235	104234	104236	
	Short	up to 17	up to 38		37 – 50	104197	104199	104198	104200	
XX-Large	Average	up to 19	up to 41	32 - 34		104217	104219	104218	104220	
	Tall	up to 22	up to 44			104237	104239	104238	104240	

#### Thigh High

sold individually Compression 20-30mmHg Top Style: Open

							Top of S	Stocking
Size	Length (cm)			Circumfer	rence (cm)	Silicone		
		(1.6.)	(10)	(	(	(cG)	Beige	Black
		(IA)	(ID)	(cB)	(cC)		REF	
	Short	up to 17	up to 71		27 – 36	41 – 55	104241	104242
Small	Average	up to 19	up to 74	18 – 20			104251	104252
	Tall	up to 22	up to 77				104261	104262
Medium	Short	up to 17	up to 71	20 – 24	29 – 39	45 – 61	104243	104244
	Average	up to 19	up to 74				104253	104254
	Tall	up to 22	up to 77				104263	104264
	Short	up to 17	up to 71		32 - 43	50 - 66	104245	104246
Large	Average	up to 19	up to 74	24 – 28			104255	104256
	Tall	up to 22	up to 77				104265	104266
	Short	up to 17	up to 71		33 – 45	56 - 71	104247	104248
X-Large	Average	up to 19	up to 74	28 – 31			104257	104258
	Tall	up to 22	up to 77				104267	104268
XX-Large	Short	up to 17	up to 71				104249	104250
	Average	up to 19	up to 74	32 - 34	37 – 50	66 - 81	104259	104260
	Tall	up to 22	up to 77				104269	104270












## ExoStrong<sup>™</sup> Measuring Guide



### **Upper Extremity**

**Circumference points:** 

cD

сС


сВ

Lengths:

- Circumference points:
- **cG** top of the arm, just below the axillary crease
- сE round the elbow joint and into the crease of the elbow
- сС narrowest point of the wrist
- сВ widest width of the palm, below the metacarpal phalangeal joints

#### Lengths:

IC-G the length between point cC and point cG







I = Length (cm) c = Circumference (cm)